



TRUEKINETIX

Packing instructions for the TrueTrainer


Step 1:

Fold in the trainer's legs and position the TrueTrainer's feet to the holes at the bottom of the box.



Step 2:

Slide the upper part with the open slot over the cassette.



Step 3:

There is a smaller box for accessories and extra protection during transport. Place it on top of the two holes and tape the box shut.

Now, your package is ready to ship.

Send the TrueTrainer to the following address:

TrueKinetix B.V.
Satellietbaan 6B
2181 MH
Hillegom
The Netherlands

